



ADDITIONAL INFORMATION

Weight	0.08 kg
Dimensions	20.5 × 8 × 3.2 cm
Carton Qty	100
Pricing Unit	Unit/1
Size	15cm x 18cm
Min Order Qty	1
ARTG ID	201258
GMDN Code	201258
Medical Device Class	Class 1s
UNSPSC Code	42311545
Tax Schedule	Taxable
Item Length	18
Item Width	15
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

FirstCare Trauma & Haemorrhage Control Bandage consolidates a sterile wound pad, robust elastic wrap, and an integral pressure applicator into a single, ready-to-deploy solution for external bleeding control. Built for pre-hospital responders, workplace first aiders, and community emergencies, it enables fast placement of direct pressure and helps maintain it through movement and handover. The all-in-one format reduces the number of separate pieces to carry and simplifies standard operating procedures when time is critical.

The sterile pad is positioned directly over the site while the wrap provides circumferential tension. A built-in pressure applicator

focuses compression exactly where it is needed to address bleeding, and a secondary dressing component helps keep the pad stable as the patient is moved. The non-adherent contact layer is designed for cleaner release at changeout, supporting low-trauma removal and efficient progression to definitive care.

To streamline application in confined spaces or stressful conditions, the system includes a closure bar that fastens without clips, pins, tape, or knots. The design also supports one-handed self-application, aiding effective use when assistance is limited. This practical construction supports consistent technique across teams, reduces training complexity, and helps responders establish targeted compression quickly with minimal equipment.

- **Direct pressure made simple:** integral applicator concentrates compression at the wound.
- **Sterile pad with secondary dressing:** helps maintain placement and pressure in transit.
- **Fast, tool-free closure:** built-in bar avoids clips, pins, tape, and knots.
- **Low-trauma removal:** non-adherent contact layer designed to lift cleanly.
- **Supports one-handed use:** aids effective application when help is limited.

The **Size 10 × 17 cm** option provides compact carry with ample pad coverage for limb and medium-area wounds, balancing portability and effective pressure application during pre-hospital response.

Finished in **Civilian White** to suit public-facing environments and quick visual checks, this option pairs the integral pressure applicator and sterile pad with a neutral colourway for versatile inclusion in general first aid kits.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Civilian Trauma & Haemorrhage Control Wound Dressing (white)