



ADDITIONAL INFORMATION

Weight	0.005 kg
Dimensions	1 × 1 × 14 cm
Pricing Unit	Each
Min Order Qty	1
ARTG ID	331117
GMDN Code	59230
Medical Device Class	Class 2a
Tax Schedule	Taxable
Item Length	3.8
Item Measurement Unit	mm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The TPAK 14 gauge x 3.25" chest decompression needle is a compact, reliable solution for treating a tension pneumothorax. The shock-resistant case ensures the needle is secure upon insertion, and the capless chamber provides instant confirmation of needle placement.