

TPAK CHEST DECOMPRESSION NEEDLE 14G X 3.25

AER-TM303



ADDITIONAL INFORMATION

Weight 0.005 kg

Dimensions $1 \times 1 \times 14 \text{ cm}$

Pricing Unit Each

Min Order Qty 1

ARTG ID 331117

GMDN Code 59230

Medical Device Class 2a

Tax Schedule Taxable

Item Length 3.8

Item Measurement Unit mm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

The TPAK 14 gauge x 3.25" chest decompression needle is a compact, reliable solution for treating a tension pneumothorax. The shock-resistant case ensures the needle is secure upon insertion, and the capless chamber provides instant confirmation of needle placement.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.