



## ADDITIONAL INFORMATION

<b>Weight</b>	0.005 kg
<b>Dimensions</b>	1 × 1 × 14 cm
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	250698
<b>GMDN Code</b>	12732
<b>Medical Device Class</b>	Class 2a
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	3.8
<b>Item Measurement Unit</b>	mm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

The TPAK 14 gauge x 3.25" chest decompression needle is a compact, reliable solution for treating a tension pneumothorax. The shock-resistant case ensures the needle is secure upon insertion, and the capless chamber provides instant confirmation of needle placement.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.