



## ADDITIONAL INFORMATION

<b>Weight</b>	0.11 kg
<b>Dimensions</b>	10.4 × 16.4 × 4 cm
<b>Carton Qty</b>	48
<b>Pricing Unit</b>	Box/100
<b>Size</b>	7.5cm x 2.5cm
<b>Min Order Qty</b>	12
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	7.5
<b>Item Width</b>	2.5
<b>Item Height</b>	7.5
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60
<b>Warranty Period (Months)</b>	0



## PRODUCT DESCRIPTION

A-CARE Detectable blue dressings provide reliable wound coverage designed for environments where contamination control matters. The bold blue colour is easy to see during visual checks, while an internal metal strip allows detection by standard metal detectors used in food processing and catering settings. This combination supports rapid identification and removal if a dressing becomes dislodged, helping reduce contamination risk across production lines and prep areas.

The washproof film backing stands up to routine handwashing and incidental moisture, maintaining a clean, low-profile seal that helps protect the site from dirt. A soft, hypo-allergenic, latex-free adhesive offers dependable hold and skin compatibility for extended wear, making these dressings suitable for sensitive skin. The flexible backing conforms to hands and fingers to limit edge-lift during movement.

Each dressing is sterile and single-use, supporting hygienic application in high-turnover workplaces. The non-stick wound pad cushions the site while minimising trauma on removal. Practical, economical and purpose-built for food handling, hospitality, and general workplace first aid, A-CARE Detectable is a dependable option for first aid kits across kitchens, production floors, and mobile service operations.

- **Strong adhesion:** secure hold that resists edge lift in busy, wet tasks.
- **Hypo-allergenic & latex-free:** suitable for sensitive skin.
- **Blue visual + metal detectable:** quick identification and detector pick-up.
- **Washproof backing:** withstands routine handwashing.
- **Sterile, single-use:** helps maintain hygiene standards in workplaces.

7.5 x 2.5 cm strips provide everyday coverage for fingers and hands, delivering a conforming fit with a low-profile seal suitable for manual tasks in food handling environments.

Packaged in **boxes of 100** to streamline replenishment for kitchens and production areas, reducing mid-shift changeovers and supporting consistent availability in first aid stations.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Standard Strip