



## ADDITIONAL INFORMATION

<b>Weight</b>	0.11 kg
<b>Dimensions</b>	4.1 × 12.8 × 10 cm
<b>Carton Qty</b>	36
<b>Pricing Unit</b>	Box/100
<b>Size</b>	7.5 x 2.5cm
<b>Min Order Qty</b>	12
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42311546
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	7.5
<b>Item Width</b>	2.5
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	9341394016224

## PRODUCT DESCRIPTION

AeroPlast™ Standard Detectable bandages are an affordable plaster, where detection is required. Its vibrant blue colour, doesn't occur naturally in food ensuring easy detection. Its unique properties also provide metal detection capabilities, where metal detectors are in use. The ultimate in high security.

### Protective Wrapper

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening. Individually wrapped for guaranteed



# AEROPLAST STANDARD DETECTABLE STRIPS 7.5 X 2.5CM BOX/100

AD3000

quality and tamper protection.

## FEATURES

- Metal detectable
- Visually detectable
- Washproof
- Conforming
- Breathable
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond™ adhesive technology
- Ideal for the manufacturing & processing sector

This product is not a substitute for direct medical attention from a healthcare practitioner.