

AEROPLAST STANDARD DETECTABLE STRIPS 7.5 X 2.5CM BOX/100

AD3000



ADDITIONAL INFORMATION

Weight	0.11 kg
Dimensions	4.1 × 12.8 × 10 cm
Carton Qty	36
Pricing Unit	Box/100
Size	7.5 x 2.5cm
Min Order Qty	12
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	7.5
ltem Width	2.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AeroPlastTM Standard Detectable bandages are an affordable plaster, where detection is required. Its vibrant blue colour, doesn't occur naturally in food ensuring easy detection. Its unique properties also provide metal detection capabilities, where metal detectors are in use. The ultimate in high security.

Protective Wrapper A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening. Individually wrapped for guaranteed quality and tamper protection.



AD3000

FEATURES

- Metal detectable
- Visually detectable
- Washproof
- Conforming
- Breathable
- Sterile
- Latex Free
- Hypoallergenic
- SupaBondTM adhesive technology
- Ideal for the manufacturing & processing sector

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.