



ADDITIONAL INFORMATION

Weight	0.06 kg
Dimensions	4.1 × 8 × 10.2 cm
Carton Qty	60
Pricing Unit	Box/50
Size	7.5 x 2.5cm
Min Order Qty	12
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC Code	42311546
Tax Schedule	Taxable
Item Length	7.5
Item Width	2.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AEROPLAST Standard Detectable plasters are designed for environments where both visual and metal detection are essential, such as food manufacturing, processing plants, commercial kitchens, and workplace canteens. The vibrant blue colour provides strong visual contrast against most foods and work surfaces, helping teams quickly identify and remove any lost plaster. A built-in metal detectable foil backing adds an extra layer of security, enabling reliable detection where metal detection equipment is in use.

Each plaster combines washproof protection with a flexible, conforming construction that moves comfortably with fingers, hands, and

other high-mobility areas. The waterproof outer layer helps protect the wound from dirt and germs while remaining air-permeable so skin can breathe. AEROPLAST Standard Detectable plasters are latex free, hypoallergenic, and designed for good skin compatibility, making them suitable for regular use and for sensitive skin in workplace and household first aid.

Under the blue detectable layer sits a super absorbent non-woven wound pad that manages exudate to help reduce skin maceration. A non-stick barrier on the pad minimises trauma on removal, supporting gentle, low-pain dressing changes. SupaBond adhesive technology delivers strong, long-lasting adhesion while remaining latex free and hypoallergenic. The adhesive is engineered to hold firmly through routine hand washing, glove changes, and active work shifts, helping reduce the need for frequent reapplication.

Protective paper and film layers support clean handling during application and help maintain the integrity of the plaster until use. The construction is designed to provide reliable wound coverage, robust adhesion strength, breathable comfort, and consistent detectability, making this plaster a dependable choice for first aid stations across food production, hospitality, healthcare support areas, and busy homes.

- **Strong adhesion:** SupaBond adhesive offers secure, long-lasting hold during active use.
- **High detectability:** vibrant blue layer and metal detectable backing support food safety systems.
- **Skin friendly:** latex free, hypoallergenic materials suitable for sensitive skin.
- **Breathable and washproof:** air-permeable, waterproof construction protects while allowing airflow.
- **Comfortable wound pad:** super absorbent, non-stick pad helps reduce maceration and supports gentle removal.

The **Size: Strips 7.5 x 2.5 cm** option provides a wide strip format that covers a larger area of fingers and hands while remaining comfortable to wear. The added width supports secure adhesion on high-movement areas such as knuckles and across joints, making it well suited to food handling tasks and other work that demands reliable coverage.

The **quantity: box of 50** offers a compact option for smaller production areas, canteens, and workplace kitchens. It supports regular replenishment without taking up excessive shelf space, helping keep blue metal detectable bandages readily available at key points of use.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.