



ADDITIONAL INFORMATION

Weight	0.06 kg
Dimensions	4.1 × 8 × 10.2 cm
Carton Qty	60
Pricing Unit	Box/50
Size	7.5 x 2.5cm
Min Order Qty	12
ARTG ID	159198
GMDN Code	44990
Medical Device Class	Class 1s
UNSPSC	42311546
Tax Schedule	Taxable
Item Length	7.5
Item Width	2.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AeroPlast™ Standard Detectable bandages are an affordable plaster, where detection is required. Its vibrant blue colour, doesn't occur naturally in food ensuring easy detection. Its unique properties also provide metal detection capabilities, where metal detectors are in use. The ultimate in high security.

Protective Wrapper

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening. Individually wrapped for guaranteed quality and tamper protection.



AEROPLAST STANDARD DETECTABLE STRIPS 7.5 X 2.5CM BOX/50

AD3050

FEATURES

- Metal detectable
- Visually detectable
- Washproof
- Conforming
- Breathable
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond™ adhesive technology
- Ideal for the manufacturing & processing sector

This product is not a substitute for direct medical attention from a healthcare practitioner.