



ADDITIONAL INFORMATION

Weight	0.02 kg
Dimensions	11.4 × 2.5 × 7.1 cm
Carton Qty	20
Pricing Unit	Each
Min Order Qty	1
ARTG ID	386216
GMDN Code	63222
Medical Device Class	Class 1
Tax Schedule	Non-Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

CPRWrap® Aid is an innovative tool that helps make CPR easier. In an emergency, just remove the product from its quick-open package, lay it on the chest of the sudden cardiac arrest (SCA) victim, and follow the simple instructions.

Made with a non-woven fabric that is eco-friendly, this is more than a CPR mask. It provides visual prompts to guide the rescuer through the AHA-recommended steps of CPR. To ensure sanitary protection, the non-woven chest overlay includes a mouth-to-mouth, one-way breathing barrier, and face shield.

- Step-by-step guide for the four AHA-recommended steps
- Visual guides ensure proper hand placement
- One-way valve mouth barrier protects from contaminants
- Non-woven (eco friendly) fabric
- Easy-to-open package for quick access in a medical emergency

- AED compatible
- Bilingual (ENG/SPAN)
- For ages 8 and up

Aid Size	Product Dimensions	Pouch Measurements	Pouch Weight
Adult	16" W x 31.5" H	3.5" W x 3.5" H	1.4 oz
Child	9.4" W x 23" H	3.5" W x 3.5" H	0.8 oz
Infant	8" W x 14.8" H	3.5" W x 3.5" H	0.4 oz

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.