

CARDIACT ALARMED OUTDOOR AED CABINET WITH LOCK 38 X 42 X 22CM

CC-60



ADDITIONAL INFORMATION

Weight 9 kg

Dimensions $47 \times 31 \times 45.5 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1

Min Order Qty 1

UNSPSC 42192416

Tax Schedule Taxable

Item Length 42

Item Width 22

Item Height 38

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

AED Alarmed Outdoor Cabinet w/lock 48cm x 47cm x 31cm Weather Rating: IP55

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.