





ADDITIONAL INFORMATION

Pricing Unit

Size 36x43x45cm **Min Order Qty** 1 **Tax Schedule** Taxable 36 **Item Length Item Width** 43 **Item Height** 45 Item Measurement Unit cm Is Serialized Item No **Is Dangerous Good** No

Each

PRODUCT DESCRIPTION

The AED Cabinet Shroud is designed to provide additional protection to AEDs and Cabinets located outdoors. Manufactured from powder-coated stainless steel the shroud provides protection for the Cabinet and AED from extreme weather conditions and provides UV protection. Keeping the AED safe and ready to use.

The shroud is suitable for our entire range of AED cabinets and our AED stand.

- Dimensions: 36x43x45 cm
- Compatibility: Designed to work with all Cabinets.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AED CABINET STAND SHROUD

CA-CS-S