



## ADDITIONAL INFORMATION

<b>Pricing Unit</b>	Each
<b>Size</b>	36x43x45cm
<b>Min Order Qty</b>	1
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	36
<b>Item Width</b>	43
<b>Item Height</b>	45
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

The AED Cabinet Shroud is designed to provide additional protection to AEDs and Cabinets located outdoors. Manufactured from powder-coated stainless steel the shroud provides protection for the Cabinet and AED from extreme weather conditions and provides UV protection. Keeping the AED safe and ready to use.

The shroud is suitable for our entire range of AED cabinets and our AED stand.

- Dimensions: 36x43x45 cm
- Compatibility: Designed to work with all Cabinets.