

AED CABINET STAND SHROUD

CA-CS-S



ADDITIONAL INFORMATION

Pricing Unit	Each
Size	36x43x45cm
Min Order Qty	1
Tax Schedule	Taxable
Item Length	36
ltem Width	43
Item Height	45
Item Measurement Unit	cm
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The AED Cabinet Shroud is designed to provide additional protection to AEDs and Cabinets located outdoors. Manufactured from powder-coated stainless steel the shroud provides protection for the Cabinet and AED from extreme weather conditions and provides UV protection. Keeping the AED safe and ready to use.

The shroud is suitable for our entire range of AED cabinets and our AED stand.

- Dimensions: 36x43x45 cm
- Compatibility: Designed to work with all Cabinets.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AED CABINET STAND SHROUD

CA-CS-S