

## CARDIACT STANDARD AED CABINET 45 X 35.5 X 14.5CM

**CC-25** 



## **ADDITIONAL INFORMATION**

Weight 3.94 kg

**Dimensions**  $48.5 \times 15.5 \times 38.6 \text{ cm}$ 

Carton Qty 1

Pricing Unit Unit/1

**Size** 450mm x 355mm x 145mm

Min Order Qty 1

UNSPSC 42192416
Tax Schedule Taxable

**Item Length** 45

Item Width 35.5

**Item Height** 14.5

Item Measurement Unit cm

**Is Serialized Item** No

Is Dangerous Good No

## **PRODUCT DESCRIPTION**

This indoor metal wall cabinet comes with Perspex front and complete with wall brackets. With its ARC approved AED signage, this cabinet ensures that your defibrillator is highly visible, increasing staff awareness and minimising lost time in time of an emergency. This attractive cabinet will also protect your AED from dust and grime.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Standard Wall Cabinet for PAD-350P/500P



## CARDIACT STANDARD AED CABINET 45 X 35.5 X 14.5CM

**CC-25**