



ADDITIONAL INFORMATION

| | |
|------------------------------|-----------------------|
| Weight | 5.76 kg |
| Dimensions | 47.9 × 20.5 × 47.2 cm |
| Carton Qty | 1 |
| Pricing Unit | Unit/1 |
| Size | 420mm x 380mm x 155mm |
| Min Order Qty | 1 |
| UNSPSC Code | 42192416 |
| Tax Schedule | Taxable |
| Item Length | 42 |
| Item Width | 38 |
| Item Height | 15.5 |
| Item Measurement Unit | cm |
| Is Serialized Item | No |
| Is Dangerous Good | No |



PRODUCT DESCRIPTION

The CC-50 is an AED wall cabinet with an alarm, strobe light, and ANZCOR approved signage. Its design allows for full visibility of the AED and is compact enough to fit in confined spaces.

- Strobe light: To attract attention to the cabinet when the door is opened
- Visible: ANZCOR approved AED signage on the sides of the cabinet
- High Quality: Premium materials with powder coated finish

- Deterrent Alarm: alerting rescuers and deterring would-be thieves
- Access: Perspex door allows full visibility of the AED
- Compact: Large enough to store most AEDs but small enough to fit in confined spaces
- Battery Powered: Uses a 9-volt battery - no hardwiring required.

Suits most conventional AED models, including HeartSine and Lifepak AEDs.
Internal Dimensions: 153 x 415 x 378mm

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

CardiAct AED Wall Cabinet with Alarm & Strobe