


CARDIACT ALARMED AED CABINET WITH STROBE LIGHT 42 X 38 X 15.5CM

CC-50





ADDITIONAL INFORMATION

Weight 5.76 kg

Dimensions $47.9 \times 20.5 \times 47.2$ cm

Carton Qty 1

Pricing Unit Unit/1

420mm x 380mm x 155mm Size

38

Min Order Qty

UNSPSC 42192416 Taxable **Tax Schedule**

Item Length 42 **Item Width**

Item Height 15.5

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

The CC-50 is an AED wall cabinet with an alarm, strobe light, and ANZCOR approved signage. Its design allows for full visibility of the AED and is compact enough to fit in confined spaces.

- Strobe light: To attract attention to the cabinet when the door is opened
- Visible: ANZCOR approved AED signage on the sides of the cabinet
- High Quality: Premium materials with powder coated finish
- Deterrent Alarm: alerting rescuers and deterring would-be thieves
- Access: Perspex door allows full visibility of the AED
- Compact: Large enough to store most AEDs but small enough to fit in confined spaces



CARDIACT ALARMED AED CABINET WITH STROBE LIGHT 42 X 38 X 15.5CM

CC-50

• Battery Powered: Uses a 9-volt battery - no hardwiring required.

Suits most conventional AED models, including HeartSine and Lifepak AEDs.

Internal Dimensions: 153 x 415 x 378mm

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

CardiAct AED Wall Cabinet with Alarm & Strobe