



ADDITIONAL INFORMATION

Weight	5.76 kg
Dimensions	47.9 × 20.5 × 47.2 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	420mm x 380mm x 155mm
Min Order Qty	1
UNSPSC	42192416
Tax Schedule	Taxable
Item Length	42
Item Width	38
Item Height	15.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The CC-50 is an AED wall cabinet with an alarm, strobe light, and ANZCOR approved signage. Its design allows for full visibility of the AED and is compact enough to fit in confined spaces.

- **Strobe light:** To attract attention to the cabinet when the door is opened
- **Visible:** ANZCOR approved AED signage on the sides of the cabinet
- **High Quality:** Premium materials with powder coated finish
- **Deterrent Alarm:** alerting rescuers and deterring would-be thieves
- **Access:** Perspex door allows full visibility of the AED
- **Compact:** Large enough to store most AEDs but small enough to fit in confined spaces



CARDIACT ALARMED AED CABINET WITH STROBE LIGHT 42 X 38 X 15.5CM

CC-50

- Battery Powered: Uses a 9-volt battery - no hardwiring required.

Suits most conventional AED models, including HeartSine and Lifepak AEDs.

CardiAct AED Wall Cabinet with Alarm & Strobe