



ADDITIONAL INFORMATION

Weight	3.06 kg
Dimensions	44.4 × 16.3 × 39.6 cm
Carton Qty	4
Pricing Unit	Unit/1
Size	External – 430mm x 380mm x 154mm, Internal – 395mm x 320mm x 145mm
Min Order Qty	1
UNSPSC	42192416
Tax Schedule	Taxable
Item Length	43
Item Width	38
Item Height	15.4
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	29341394027217

PRODUCT DESCRIPTION

This sturdy and rugged water and airtight case is ideal for exceptionally humid or rugged conditions. This customisable case is essential for all marine and outdoor applications.

Rated to IP67 this case can also be equipped with the Desiccant reusable Silica Gel to remove any moisture that may be present in the atmosphere after opening the case.

Fits all **LifePak CR2 Defibrillators** (CR-2, CR-2-A, CR-2E, CR-2E-A).

Large Waterproof Tough AED Case for all LifePak CR2 Defibrillators (CR-2, CR-2-A, CR-2E, CR-2E-A).



CARDIACT LARGE WATERPROOF TOUGH AED CASE 43 X 38 X 15.4CM (YELLOW)

CW-20