



## ADDITIONAL INFORMATION

<b>Weight</b>	3.06 kg
<b>Dimensions</b>	44.4 × 16.3 × 39.6 cm
<b>Carton Qty</b>	4
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	External – 430mm x 380mm x 154mm, Internal – 395mm x 320mm x 145mm
<b>Min Order Qty</b>	1
<b>UNSPSC Code</b>	42192416
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	43
<b>Item Width</b>	38
<b>Item Height</b>	15.4
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

This sturdy and rugged water and airtight case is ideal for exceptionally humid or rugged conditions. This customisable case is essential for all marine and outdoor applications. Rated to IP67 this case can also be equipped with the Desiccant reusable Silica Gel to remove any moisture that may be present in the atmosphere after opening the case. Fits all **LifePak CR2 Defibrillators** (CR-2, CR-2-A, CR-2E, CR-2E-A).

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Large Waterproof Tough AED Case for all LifePak CR2 Defibrillators (CR-2, CR-2-A, CR-2E, CR-2E-A).



# CARDIAC LARGE WATERPROOF TOUGH AED CASE 43 X 38 X 15.4CM (YELLOW)

CW-20