



ADDITIONAL INFORMATION

| | |
|------------------------------|------------------|
| Weight | 0.4 kg |
| Dimensions | 24 × 2.5 × 32 cm |
| Carton Qty | 12 |
| Pricing Unit | Wrap/12 |
| Size | 10cm x 4.5M |
| Min Order Qty | 1 |
| ARTG ID | 236794 |
| GMDN Code | 34864 |
| Medical Device Class | Class 1 |
| UNSPSC Code | 42311500 |
| Tax Schedule | Taxable |
| Item Length | 450 |
| Item Width | 10 |
| Item Measurement Unit | cm |
| Is Medical Device | Yes |
| Is Serialized Item | No |
| Is Dangerous Good | No |
| Shelf Life (Months) | 60 |



PRODUCT DESCRIPTION

AeroBan Cohesive Bandage Blue is designed to provide comfortable, dependable support without sticking to skin, hair, or clothing. The cohesive material clings to itself, creating a secure wrap that stays in place once applied while remaining easy to remove. This makes it an excellent choice for securing dressings, supporting joints and soft-tissue injuries, or providing light compression in everyday first aid situations.

The blue colour offers high visibility, which is especially useful in food handling, catering, and workplace environments where easy identification of lost dressings is important. The bandage is non-sterile and latex free, supporting use on a wide range of users, including those who are sensitive to natural rubber latex. Breathable construction helps maintain airflow to the covered area, promoting comfort during extended wear.

AeroBan Cohesive Bandage Blue is engineered to deliver compressive support while remaining flexible and conforming around contours such as wrists, ankles, fingers, and limbs. Because it holds its position without the need for clips or fasteners, it simplifies application and helps reduce the risk of the bandage unravelling during normal movement. The material can be torn by hand, allowing quick adjustment to the required length without scissors, which is ideal for fast response in first aid kits, workplaces, sports settings, and home care.

With its combination of premium quality construction, cohesive self-adhesion, and latex-free materials, AeroBan Cohesive Bandage Blue is a practical addition to first aid kits, occupational health supplies, and everyday medical cupboards. It is suitable for retaining dressings, providing gentle compression, and offering light support to minor strains or sprains, helping users stay mobile and comfortable while protected.

- **Cohesive design:** sticks to itself, not to skin or hair, for comfortable wear and easy removal.
- **Supportive and compressive:** helps secure dressings and provide light support to joints and soft tissues.
- **Latex free and breathable:** suitable for users with latex sensitivity and designed for comfort during extended wear.
- **Hand-tearable material:** can be torn to length without scissors for quick, convenient application.
- **High-visibility blue:** easy to see in workplaces and food preparation environments.

The **10 cm x 4.5 m** size provides a wide bandage with generous wrap length, ideal for covering larger areas of the lower leg, knee, or forearm. Its broader profile helps achieve stable, even compression and dressing retention where extra coverage and support are required.

Supplied as a single 10 cm x 4.5 m blue cohesive wrap, this option delivers generous coverage for larger areas such as the lower leg, knee, or forearm. The wider width helps achieve even compression and secure dressing retention while still fitting neatly into workplace and home first aid kits.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.