



ADDITIONAL INFORMATION

Weight	0.12 kg
Dimensions	13 × 5.5 × 19 cm
Carton Qty	48
Pricing Unit	Wrap/12
Size	2.5cm x 4.5M
Min Order Qty	1
ARTG ID	236794
GMDN Code	34864
Medical Device Class	Class 1
UNSPSC Code	42311500
Tax Schedule	Taxable
Item Length	450
Item Width	2.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60



PRODUCT DESCRIPTION

AeroBan Cohesive Bandage Blue is designed to provide comfortable, dependable support without sticking to skin, hair, or clothing. The cohesive material clings to itself, creating a secure wrap that stays in place once applied while remaining easy to remove. This makes it an excellent choice for securing dressings, supporting joints and soft-tissue injuries, or providing light compression in everyday first aid situations.

The blue colour offers high visibility, which is especially useful in food handling, catering, and workplace environments where easy identification of lost dressings is important. The bandage is non-sterile and latex free, supporting use on a wide range of users, including those who are sensitive to natural rubber latex. Breathable construction helps maintain airflow to the covered area, promoting comfort during extended wear.

AeroBan Cohesive Bandage Blue is engineered to deliver compressive support while remaining flexible and conforming around contours such as wrists, ankles, fingers, and limbs. Because it holds its position without the need for clips or fasteners, it simplifies application and helps reduce the risk of the bandage unravelling during normal movement. The material can be torn by hand, allowing quick adjustment to the required length without scissors, which is ideal for fast response in first aid kits, workplaces, sports settings, and home care.

With its combination of premium quality construction, cohesive self-adhesion, and latex-free materials, AeroBan Cohesive Bandage Blue is a practical addition to first aid kits, occupational health supplies, and everyday medical cupboards. It is suitable for retaining dressings, providing gentle compression, and offering light support to minor strains or sprains, helping users stay mobile and comfortable while protected.

- **Cohesive design:** sticks to itself, not to skin or hair, for comfortable wear and easy removal.
- **Supportive and compressive:** helps secure dressings and provide light support to joints and soft tissues.
- **Latex free and breathable:** suitable for users with latex sensitivity and designed for comfort during extended wear.
- **Hand-tearable material:** can be torn to length without scissors for quick, convenient application.
- **High-visibility blue:** easy to see in workplaces and food preparation environments.

The **Size 2.5 cm x 4.5 m** option offers a slim, versatile wrap that is ideal for fingers, wrists, paws, and smaller limbs, providing neat coverage and light compression while remaining easy to carry in compact first aid kits and bags.

Supplied as a cohesive bandage measuring **2.5 cm x 4.5 m (stretched)**, this option balances compact width with ample length for multiple turns around joints or limbs, making it suitable for single applications and convenient replenishment in household or workplace first aid kits.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Cohesive Bandage 2.5cm