

AEROBAN BLUE COHESIVE BANDAGE 5.0CM X 4.5M WRAP/12

ABCH50





ADDITIONAL INFORMATION

Weight 0.2 kg

Dimensions 17 × 5.5 × 2.5 cm

Carton Qty 24

Pricing Unit Wrap/12

Size 5cm x 4.5M

Min Order Qty 1

ARTG ID 236794 **GMDN** Code 34864

Medical Device Class Class 1

UNSPSC 42311500

Tax Schedule Taxable

Item Length 450

Item Width 5

Item Measurement Unit cm

Is Medical Device Yes Is Serialized Item

No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

AEROBAN Cohesive Bandages are a premium quality bandage that has the unique ability to stick to itself but not to skin

Cohesive

• This bandage provides support by adhering to itself, and not to skin or hair. It also will hold itself in position without the use of bandage clips



AEROBAN BLUE COHESIVE BANDAGE 5.0CM X 4.5M WRAP/12

ABCH50

Easily adjusted to size

• AEROBAN Cohesive Bandages can be easily be torn to size, eliminating the need to use scissors

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Cohesive Bandage 5cm