



ADDITIONAL INFORMATION

Weight	0.33 kg
Dimensions	30 × 15 × 5 cm
Carton Qty	16
Pricing Unit	Wrap/12
Size	7.5cm x 4.5M
Min Order Qty	1
ARTG ID	236794
GMDN Code	34864
Medical Device Class	Class 1
UNSPSC	42311500
Tax Schedule	Taxable
Item Length	450
Item Width	7.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394001374

PRODUCT DESCRIPTION

AEROBAN Cohesive Bandages are a premium quality bandage that has the unique ability to stick to itself but not to skin

Cohesive

- This bandage provides support by adhering to itself, and not to skin or hair. It also will hold itself in position without the use of bandage clips



AEROBAN BLUE COHESIVE BANDAGE 7.5CM X 4.5M WRAP/12

ABCH75

Easily adjusted to size

- AEROBAN Cohesive Bandages can be easily be torn to size, eliminating the need to use scissors

Cohesive Bandage 7.5cm