

## AEROBAN TAN COHESIVE BANDAGE 10CM X 4.5M WRAP/12

**ABCT100** 



#### **ADDITIONAL INFORMATION**

Weight 0.4 kg

**Dimensions**  $5.5 \times 24 \times 32 \text{ cm}$ 

Carton Qty 12

Pricing Unit Wrap/12

**Size** 10cm x 4.5M

Min Order Qty 1

ARTG ID 236794
GMDN Code 34864
Medical Device Class 1

**UNSPSC** 42311500

Tax Schedule Taxable

Item Length 450

Item Width 10

Item Measurement Unit cm

Is Medical Device Yes

**Is Serialized Item** No

Is Dangerous Good No

Shelf Life (Months) 60

#### **PRODUCT DESCRIPTION**

AEROBAN Tan Cohesive Bandages are a premium quality bandage with the unique ability to stick to itself but not to skin. Available in a tan colour for discrete application.

#### Cohesive

• This bandage provides support by adhering to itself, and not to skin or hair. It also will hold itself in position without the use of bandage clip

### Easily adjusted to Size



# AEROBAN TAN COHESIVE BANDAGE 10CM X 4.5M WRAP/12

**ABCT100** 

AEROBAN Cohesive Bandages can be easily be torn to size, eliminating the need to use scissors.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Tan Cohesive Bandage 10cm