



ADDITIONAL INFORMATION

| | |
|------------------------------|------------------|
| Weight | 0.2 kg |
| Dimensions | 5.5 × 17 × 25 cm |
| Carton Qty | 24 |
| Pricing Unit | Wrap/12 |
| Size | 5cm x 4.5M |
| Min Order Qty | 1 |
| ARTG ID | 236794 |
| GMDN Code | 34864 |
| Medical Device Class | Class 1 |
| UNSPSC | 34864 |
| Tax Schedule | Taxable |
| Item Length | 450 |
| Item Width | 5 |
| Item Measurement Unit | cm |
| Is Serialized Item | No |
| Is Dangerous Good | No |
| Shelf Life (Months) | 60 |

PRODUCT DESCRIPTION

AEROBAN Tan Cohesive Bandages are a premium quality bandage with the unique ability to stick to itself but not to skin. Available in a tan colour for discrete application.

Cohesive

- This bandage provides support by adhering to itself, and not to skin or hair. It also will hold itself in position without the use of bandage clip

Easily adjusted to Size

AEROBAN Cohesive Bandages can be easily be torn to size, eliminating the need to use scissors.



AEROBAN TAN COHESIVE BANDAGE 5.0CM X 4.5M WRAP/12

ABCT50

Tan Cohesive Bandage 5cm