

AEROBAN TAN COHESIVE BANDAGE 7.5CM X 4.5M WRAP/12

ABCT75



ADDITIONAL INFORMATION

Weight 0.33 kg

Dimensions $5.5 \times 23 \times 26 \text{ cm}$

Carton Qty 16

Pricing Unit Wrap/12

Min Order Qty 1

ARTG ID 236794
GMDN Code 34864
Medical Device Class 1

UNSPSC 42311500

Tax Schedule Taxable

Is Medical DeviceYesIs Serialized ItemNoIs Dangerous GoodNoShelf Life (Months)60

PRODUCT DESCRIPTION

AEROBAN Tan Cohesive Bandages are a premium quality bandage with the unique ability to stick to itself but not to skin. Available in a tan colour for discrete application.

Cohesive

• This bandage provides support by adhering to itself, and not to skin or hair. It also will hold itself in position without the use of bandage clip

Easily adjusted to Size

AEROBAN Cohesive Bandages can be easily be torn to size, eliminating the need to use scissors.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AEROBAN TAN COHESIVE BANDAGE 7.5CM X 4.5M WRAP/12

ABCT75

Tan Cohesive Bandage 7.5cm