



ADDITIONAL INFORMATION

Weight	0.006 kg
Dimensions	2.5 × 3.5 × 8.5 cm
Pricing Unit	Unit/1
Size	Size 01
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AEROFORM Tubular Bandage Finger Applicator is a plastic applicator used to hold the AB0120 finger tubular bandage in place, which are suited to treatment of sports injuries, dislocations and sprains.

Plastic Applicator