



## ADDITIONAL INFORMATION

<b>Weight</b>	0.037 kg
<b>Dimensions</b>	4 × 7.5 × 15.5 cm
<b>Carton Qty</b>	48
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	Size D x 1M
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	362410
<b>GMDN Code</b>	61652
<b>Medical Device Class</b>	Class 1
<b>UNSPSC Code</b>	42311506
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	100
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No



## PRODUCT DESCRIPTION

AEROFORM Elastic Tubular Bandage delivers flexible compression and soft-tissue support in a seamless, easy-to-apply tube. The elasticised knit conforms smoothly to limbs and joints without the need for clips or tapes, providing even, circumferential pressure that helps reduce swelling and protects dressings during movement. It can be trimmed to the required length and layered to adjust firmness, making it a practical choice for clinical environments, workplaces, sports settings, and home first aid.

The breathable, latex-free construction promotes skin comfort during extended wear while maintaining consistent support through activity. A smooth, low-friction surface slides easily beneath clothing or protective gear to minimise catching and disturbance, helping applications stay neat and secure during daily tasks.

Common applications include securing primary dressings, providing even compression during rehabilitation, and offering extra reassurance for weak or aching joints. The tubular design enables rapid application and a tidy, professional finish when time and consistency matter.

- **Elasticised tubular design:** delivers even, circumferential compression for secure support.
- **Cut-to-length convenience:** trim and layer to tailor firmness and coverage.
- **Breathable, latex-free knit:** promotes comfort while you move.
- **Non-sterile format:** ideal for routine support and dressing retention.

The **Size D Large Limbs (7.5 cm width)** option is suited to larger forearms and calves, providing broader coverage while still conforming closely to the limb. This size offers comfortable, even compression across a wider area, helping to keep dressings tidy and deliver reassuring support during everyday movement or light activity.

This option is supplied as a **1 metre length**, allowing the bandage to be cut quickly to the exact coverage needed for each application. The generous length supports easy customisation and folding back over the limb for firmer compression, helping to keep dressings neat and provide reliable support for larger forearms and calves.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Size B Tubular Bandage 10M