

ABD10





## **ADDITIONAL INFORMATION**

Weight 0.065 kg

**Dimensions**  $13.5 \times 13.5 \times 0.5 \text{ cm}$ 

Carton Qty 144

Pricing Unit Unit/1

Size 10cm x 10cm

Min Order Qty 36

**ARTG ID** 221543 **GMDN Code** 47764

Medical Device Class 2a

**UNSPSC** 42311513

**Tax Schedule** Non-Taxable

Item Length10Item Width10Item Height13.5

Item Measurement Unit cm

Is Medical Device Yes
Is Serialized Item No

**Is Dangerous Good** No

Shelf Life (Months) 60

## PRODUCT DESCRIPTION

AEROBURN Burn Gel Dressings offer immediate and rapid cooling helping to prevent further tissue damage, relieve pain and aid in the healing process.

## **Hydrates the Burn Site**

• AeroBurn Burn Gel Dressings hydrate the burn site giving relief from pain and encouraging faster healing.



## **AEROBURN BURN GEL DRESSING 10 X 10CM**

ABD10

Non-adherent, Non-to	adhe	rent.	Non	-toxic
----------------------	------	-------	-----	--------

Non-adherent, Non-toxic
<ul> <li>AeroBurn Burn Gel Dressings are non adherent and non-toxic, ensuring minimal irritation and maximum relief for the wound.</li> </ul>
• AeroBurn Burn Gel Dressings offer more effective relief for burn wounds than traditional methods, as well as being more cost effective. The burn dressing can simply remain in situ for up to 2 hours, providing relief with minimal attention or redressing or cooling using running water.
Range of Sizes  • AeroBurn Burn Gel Dressings are available in sizes 10cm x 10cm, 20cm x 20cm and a face dressing 30cm x 40cm and larger 40cm x 60cm dressing.
This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.
Burn Gel Infused Dressing