



ABD10





ADDITIONAL INFORMATION

Weight 0.065 kg

Dimensions $13.5 \times 13.5 \times 0.5 \text{ cm}$

Carton Qty 144

Pricing Unit Unit/1

Size 10cm x 10cm

Min Order Qty 36

ARTG ID 221543 **GMDN Code** 47764

Medical Device Class 2a

UNSPSC 42311513

Tax Schedule Non-Taxable

Item Length10Item Width10Item Height13.5

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

AeroBurn Burn Gel Dressing is a sterile, gel-impregnated pad purpose-built for first aid treatment of minor burns and scalds in workplaces, community settings, vehicles and home kits. Its water-based hydrogel helps cool the burn surface and soothe discomfort on contact, while the soft non-woven carrier lays smoothly over curved areas for quick, confident coverage. The single-use format supports best practice first response and straightforward disposal after reassessment.

Designed as a protective covering, the dressing functions as a mechanical barrier during the critical early phase of care. The gel-infused construction helps manage minimal to moderate exudate, and the low-adherent interface supports gentle removal when it's time to review the injury. Because the hydrogel is water soluble, it can be rinsed away easily with clean water to support follow-up care and redressing.



AEROBURN BURN GEL DRESSING 10 X 10CM

ABD10

Each dressing is sterile and sealed in an easy-open pouch to aid fast access when seconds matter. Trusted by organisations that prioritise readiness and reliability, AeroBurn Burn Gel Dressing integrates seamlessly into first aid stations and mobile kits across kitchens, workshops, warehouses, classrooms and field vehicles.

- Cools and soothes: water-based hydrogel helps cool burn surfaces and ease discomfort
- Protective barrier: sterile covering acts as a mechanical barrier during initial care
- **Soft non-woven carrier:** conforms comfortably for rapid, even coverage
- Low-adherent contact: designed for gentle removal to support reassessment
- Water-soluble gel: rinses off with clean water to simplify follow-up care

Size – 10 × 10 cm: standard size, combining a soft non-woven carrier with cooling hydrogel to drape smoothly and support even contact during initial burn care.

Supplied as a **single sterile pouch** containing one **10 cm** × **10 cm** gel-infused dressing. Compact and kit-ready for quick cooling coverage of minor burns before follow-up care.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Burn Gel Infused Dressing