



## ADDITIONAL INFORMATION

<b>Weight</b>	0.065 kg
<b>Dimensions</b>	13.5 × 13.5 × 0.5 cm
<b>Carton Qty</b>	144
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	10cm x 10cm
<b>Min Order Qty</b>	36
<b>ARTG ID</b>	221543
<b>GMDN Code</b>	47764
<b>Medical Device Class</b>	Class 2a
<b>UNSPSC</b>	42311513
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	10
<b>Item Width</b>	10
<b>Item Height</b>	13.5
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60

## PRODUCT DESCRIPTION

AEROBURN Burn Gel Dressings offer immediate and rapid cooling helping to prevent further tissue damage, relieve pain and aid in the healing process.

### Hydrates the Burn Site

- AeroBurn Burn Gel Dressings hydrate the burn site giving relief from pain and encouraging faster healing.

## **Non-adherent, Non-toxic**

- AeroBurn Burn Gel Dressings are non adherent and non-toxic, ensuring minimal irritation and maximum relief for the wound.

## **Easy to Use**

- AeroBurn Burn Gel Dressings offer more effective relief for burn wounds than traditional methods, as well as being more cost effective. The burn dressing can simply remain in situ for up to 2 hours, providing relief with minimal attention or redressing or cooling using running water.

## **Range of Sizes**

- AeroBurn Burn Gel Dressings are available in sizes 10cm x 10cm, 20cm x 20cm and a face dressing 30cm x 40cm and larger 40cm x 60cm dressing.

Burn Gel Infused Dressing