



## ADDITIONAL INFORMATION

<b>Weight</b>	0.17 kg
<b>Dimensions</b>	27 × 23 × 0.5 cm
<b>Carton Qty</b>	48
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	75cm x 220cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	218614
<b>GMDN Code</b>	36093
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42311545
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	220
<b>Item Width</b>	75
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	9341394004078

## PRODUCT DESCRIPTION

AEROBURN Burn Sheets are a sterile sheet available in multiple sizes for covering and protecting burn wounds.

### Provides user comfort and wound protection

- The AEROBURN Low Adherent Burn Sheet is coated on one side with a low-adherent coating which minimises trauma on removal. Polyethylene Burn Sheets form a protective and non-adherent shield over burn wounds to prevent infection and protect the wound site.

## Premium Quality Materials

- AEROBURN burn sheets use premium quality materials which provide a secure barrier to infection of the wound site. AEROBURN burn sheets conforms to body contours for increased user convenience and patient comfort.

## Latex Free

- AEROBURN burn sheets are latex free, minimising risk of patient reaction or discomfort.

## Sterile

- AEROBURN burn sheets are sterile, minimising risk of secondary infection

## Convenient Application

- AEROBURN burn sheets can be easily folded or cut to size for application on all sizes of wound sites.

Low Adherent Burn Sheet