



ADDITIONAL INFORMATION

Weight	0.06 kg
Dimensions	26 × 17 × 0.5 cm
Carton Qty	96
Pricing Unit	Unit/1
Size	75cm x 75cm
Min Order Qty	1
ARTG ID	218614
GMDN Code	36093
Medical Device Class	Class 1s
UNSPSC	42311545
Tax Schedule	Taxable
Item Length	75
Item Width	75
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AEROBURN Burn Sheets are a sterile sheet available in multiple sizes for covering and protecting burn wounds.

Provides user comfort and wound protection

- The AEROBURN Low Adherent Burn Sheet is coated on one side with a low-adherent coating which minimises trauma on removal. Polyethylene Burn Sheets form a protective and non-adherent shield over burn wounds to prevent infection and protect the wound site.

Premium Quality Materials

- AEROBURN burn sheets use premium quality materials which provide a secure barrier to infection of the wound site. AEROBURN burn sheets conforms to body contours for increased user convenience and patient comfort.

Latex Free

- AEROBURN burn sheets are latex free, minimising risk of patient reaction or discomfort.

Sterile

- AEROBURN burn sheets are sterile, minimising risk of secondary infection

Convenient Application

- AEROBURN burn sheets can be easily folded or cut to size for application on all sizes of wound sites.

Low Adherent Burn Sheet