

# AEROBURN POLYETHYLENE BURN SHEET 20 X 20CM BOX/25

**ABS2020** 





Weight 0.1 kg

**Dimensions**  $14.5 \times 8.5 \times 6 \text{ cm}$ 

Carton Qty 24

Pricing Unit Box/25

**Size** 20cm x 20cm

Min Order Qty 1

**ARTG ID** 218614 **GMDN Code** 36093

Medical Device Class 1s

**UNSPSC** 42311545

Tax Schedule Taxable

Item Length 20

Item Width 20

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

Shelf Life (Months) 60

**Is Dangerous Good** 





#### **PRODUCT DESCRIPTION**

#### Provides user comfort and wound protection

• The AEROBURN Low Adherent Burn Sheet is coated on one side with a low-adherent coating which minimises trauma on removal. Polyethylene Burn Sheets form a protective and non-adherent shield over burn wounds to prevent infection and protect the wound site.

No



## AEROBURN POLYETHYLENE BURN SHEET 20 X 20CM BOX/25

**ABS2020** 

Premium Quality Materials		
	<ul> <li>AEROBURN burn sheets use premium quality materials which provide a secure barrier to infection of the wound site.</li> <li>AEROBURN burn sheets conforms to body contours for increased user convenience and patient comfort.</li> </ul>	
Latex Free	<ul> <li>AEROBURN burn sheets are latex free, minimising risk of patient reaction or discomfort.</li> </ul>	
Sterile	AEROBURN burn sheets are sterile, minimising risk of secondary infection	



#### **AEROBURN POLYETHYLENE BURN SHEET 20 X 20CM BOX/25**

020

	ABS20
Convenient Application	
<ul> <li>AEROBURN burn sheets can be easily folded or cut to size for application on all sizes of wound sites.</li> </ul>	
Range of Sizes	
Latex-free square and rectangular polyethylene burn sheets, available in the following sizes:	
○ 10cm x 10cm	
○ 20cm x 20cm	
∘ 60cm x 90cm	
o docin x 90cm	

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



## AEROBURN POLYETHYLENE BURN SHEET 20 X 20CM BOX/25

**ABS2020** 

Polyethylene Burn Sheet