

AEROBURN POLYETHYLENE BURN SHEET 60 X 90CM BOX/25

ABS6090





ADDITIONAL INFORMATION

Weight 0.7 kg

Dimensions $12 \times 21 \times 12 \text{ cm}$

Carton Qty 12

Pricing Unit Box/25

Size 60cm x 90cm

Min Order Qty 1

ARTG ID 218614 **GMDN Code** 36093

Medical Device Class 1s

UNSPSC 42311545

Tax Schedule Taxable

Item Length 90

Item Width 60

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

Provides user comfort and wound protection

• The AEROBURN Low Adherent Burn Sheet is coated on one side with a low-adherent coating which minimises trauma on removal. Polyethylene Burn Sheets form a protective and non-adherent shield over burn wounds to prevent infection and protect the wound site.



AEROBURN POLYETHYLENE BURN SHEET 60 X 90CM BOX/25

ABS6090

Premium	Quality	v Mate	rials
I I CIIII GII	. Quuiit	, iviace	ı ıuıs

Premium	Quality Materials
	 AEROBURN burn sheets use premium quality materials which provide a secure barrier to infection of the wound site AEROBURN burn sheets conforms to body contours for increased user convenience and patient comfort.
Latex Fre	e • AEROBURN burn sheets are latex free, minimising risk of patient reaction or discomfort.
Sterile	AEROBURN burn sheets are sterile, minimising risk of secondary infection
Convenie	nt Application • AEROBURN burn sheets can be easily folded or cut to size for application on all sizes of wound sites.

Range of Sizes

Latex-free square and rectangular polyethylene burn sheets, available in the following sizes:

- o 10cm x 10cm
- o 20cm x 20cm
- o 60cm x 90cm



AEROBURN POLYETHYLENE BURN SHEET 60 X 90CM BOX/25

ABS6090

This p	product is not a substitute for d	irect medical attention	from a healthcare	practitioner. Alway	s read the label a	and follow the
direct	ions for use.					

Polyethylene Burn Sheet