



## ADDITIONAL INFORMATION

<b>Weight</b>	0.7 kg
<b>Dimensions</b>	12 × 21 × 12 cm
<b>Carton Qty</b>	12
<b>Pricing Unit</b>	Box/25
<b>Size</b>	60cm x 90cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	218614
<b>GMDN Code</b>	36093
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311545
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	90
<b>Item Width</b>	60
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	60



## PRODUCT DESCRIPTION

AeroBurn AEROBURN Polyethylene Burn Sheet is a sterile, low-adherent trauma covering designed for rapid protection of burn injuries in first aid and clinical settings. As part of the Burn Treatments: Dressing category, these polyethylene (PE) sheets create a non-adherent barrier over the wound to help shield the site from external contamination while minimising trauma during removal. Clear material allows ongoing visual inspection of the affected area, supporting appropriate care decisions during transport and



# AEROBURN POLYETHYLENE BURN SHEET 60 X 90CM BOX/25

ABS6090

handover.

**Low-Adherent Protection** Each sheet is coated on one side with a low-adherent finish that reduces sticking to delicate, compromised skin. This helps minimise pain and tissue disruption at dressing changes while maintaining a clean interface over the wound surface.

**Sterile and Single Use** Supplied sterile (ethylene oxide sterilisation) and intended for single use, the sheets support infection prevention protocols in pre-hospital, workplace, and clinical environments. Tamper-evident pouch packaging with clear “Peel to Open” instructions promotes fast, aseptic access when seconds matter.

**Latex Free and Patient Friendly** Construction is latex free to reduce the risk of latex sensitivity. Soft, flexible PE conforms comfortably to body contours on hands, arms, torso, and other curved areas, helping maintain an even protective layer during movement.

**Cut-to-Size Convenience** The sheets can be folded or trimmed to suit different wound geometries and coverage needs. This flexibility helps responders adapt quickly to the injury size while conserving supplies and minimising waste.

**Quality Assured** Product is sized for efficient first aid kit replenishment and bulk procurement. This reliable range equips responders with a consistent, easy-to-use barrier dressing for emergency burn protection from the scene of injury through to clinical assessment. It is suitable for inclusion in workplace kits, emergency vehicles, and community first aid programs where dependable, sterile burn coverage is required.

Supplied as a **box of 25 sheets measuring 60 × 90 cm**, ideal for bulk restocking of workplace kits and emergency vehicles, reducing changeovers and ensuring reliable coverage for larger areas.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Polyethylene Burn Sheet