

AEROCREPE ELASTIC CREPE BANDAGE 10CM X 4M WRAP/12

ACE100





ADDITIONAL INFORMATION

Weight 0.036667 kg

Dimensions $20 \times 27 \times 4.5 \text{ cm}$

Carton Qty 20

Pricing Unit Wrap/12

Size 10cm x 4M

Min Order Qty 1

ARTG ID 236934 **GMDN Code** 10279

Medical Device Class 1

UNSPSC 42311500

Tax Schedule Taxable

Item Length 400

Item Width 10

Item Measurement Unit cm

Is Medical Device Yes

Is Serialized Item

Is Dangerous Good No

No

Shelf Life (Months) 60

PRODUCT DESCRIPTION

AeroCrepe elastic crepe bandages provide dependable support, retention and light compression for everyday first aid and clinical settings. Built from cotton crepe with elastic yarns, the bandage stretches comfortably to conform around limbs and joints while allowing normal movement. This balance of elasticity, softness and breathability makes it suitable for holding wound dressings in place, providing localised pressure, or giving supportive wrap to minor strains and swelling during short-term recovery.

The construction is non-sterile and latex free, supporting broad workplace and home use. Fray-resistant edges help maintain integrity during application and wear, while a slip-resistant texture assists with secure layering as the bandage is wrapped.

Designed for reliability in first aid kits, workplace cabinets and healthcare facilities, AeroCrepe presents a professional finish that is easy to apply and comfortable for the wearer. The material's air-permeable structure helps skin remain comfortable over extended periods, and the elastic return gives consistent compression without excessive bulk. From routine dressing retention to light



AEROCREPE ELASTIC CREPE BANDAGE 10CM X 4M WRAP/12

ACE100

compression techniques, the bandage helps users respond quickly and confidently to common injuries.

- Elastic comfort and support: stretches to conform and provide controlled, local pressure
- Breathable cotton crepe: helps keep skin comfortable during wear
- Fray-resistant edges: maintains bandage integrity during application
- Slip-resistant texture: assists in secure wrapping and layer stability
- Latex free and non-sterile: suitable for general first aid use
- Clear storage guidance: symbols indicate keep dry, avoid sunlight, and temperature limits

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Elastic Crepe Bandage 10cm