

AEROCREPE ELASTIC CREPE BANDAGE 2.5CM X 4M WRAP/12

ACE25





Weight	0.13 kg
Dimensions	27 × 5 × 4.5 cm
Carton Qty	60
Pricing Unit	Wrap/12
Size	2.5cm x 4M
Min Order Qty	1
ARTG ID	236934
GMDN Code	10279
Medical Device Class	Class 1
UNSPSC	42311500
Tax Schedule	Taxable
Item Length	400
ltem Width	2.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

ADDITIONAL INFORMATION

PRODUCT DESCRIPTION

The AEROCREPE[™] Crepe Bandage is a high-quality cotton crepe-style bandage typically recommended for light to medium injury support requirements.

The crepe bandage comes in three different types (light, medium and elastic) to cater to the retention and protection of wound dressings, as well as provide assistance in standard immobilisation and compression bandage techniques.

Perform immobilisation and compression techniques with elastic bandages

The AEROCREPE[™] Elastic Cotton Crepe Bandage has a durable, elastic conforming nature, making it suitable for injury immobilisation and compression bandage techniques.

This elastic bandage can serve to compress and add pressure to the area around an open wound to help stem blood flow, especially

while awaiting treatment. The *AEROCREPETM* elastic bandage can also be used to immobilise the area around an injury, especially if the injury requires additional time and/or a lack of use to fully recover.

Comprehensive size ranges for all needs



AEROCREPE ELASTIC CREPE BANDAGE 2.5CM X 4M WRAP/12

ACE25

Each AEROCREPE[™] cotton bandage roll comes in a length of 4 metres, with a variety of widths to suit all wound and injury requirements:

AEROCREPE[™] Light Cotton Crepe Bandage 5cm, 7.5cm, 10cm, 15cm

AEROCREPE[™] Medium Cotton Crepe Bandage 2.5cm, 5cm, 7.5cm, 10cm, 15cm

AEROCREPE[™] Elastic Cotton Crepe Bandage 2.5cm, 5cm, 7.5cm, 10cm, 15cm

It's recommended to have a selection of AEROCREPE[™] cotton bandage rolls to be prepared to address any relevant injuries. **Easy to open wrapper**

All AEROCREPE[™] cotton crepe bandages have an easy-open tab built into the wrapper of the bandage to allow for quick and easy use during an emergency, while keeping the contents safe and clean when unopened.

Contraindications and Precautions

Any wound and gauze dressings underneath cotton compression bandages should be changed daily, or more often, according to the condition of the wound.

Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This Crepe Bandage product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Elastic Crepe Bandage 2.5cm