

# AEROCREPE ELASTIC CREPE BANDAGE 5CM X 4M WRAP/12

ACE50





### **ADDITIONAL INFORMATION**

Weight 0.24 kg

**Dimensions**  $27 \times 10 \times 4.5 \text{ cm}$ 

Carton Qty 40

Pricing Unit Wrap/12

Size 5cm x 4M

Min Order Qty 1

**ARTG ID** 236934 **GMDN Code** 10279

Medical Device Class 1

**UNSPSC** 42311500

**Tax Schedule** Taxable

Item Length 400

**Item Width** 5

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

**Is Dangerous Good** No

Shelf Life (Months) 60

### **PRODUCT DESCRIPTION**

The AEROCREPE<sup>™</sup> Crepe Bandage is a high-quality cotton crepe-style bandage typically recommended for light to medium injury support requirements.

The crepe bandage comes in three different types (light, medium and elastic) to cater to the retention and protection of wound dressings, as well as provide assistance in standard immobilisation and compression bandage techniques.

# Perform immobilisation and compression techniques with elastic bandages

The  $AEROCREPE^{TM}$  Elastic Cotton Crepe Bandage has a durable, elastic conforming nature, making it suitable for injury immobilisation and compression bandage techniques.



# AEROCREPE ELASTIC CREPE BANDAGE 5CM X 4M WRAP/12

ACE50

This elastic bandage can serve to compress and add pressure to the area around an open wound to help stem blood flow, especially while awaiting treatment. The  $AEROCREPE^{TM}$  elastic bandage can also be used to immobilise the area around an injury, especially if the injury requires additional time and/or a lack of use to fully recover.

## Comprehensive size ranges for all needs

Each AEROCREPE<sup>™</sup> cotton bandage roll comes in a length of 4 metres, with a variety of widths to suit all wound and injury requirements:

AEROCREPE<sup>™</sup> Light Cotton Crepe Bandage 5cm, 7.5cm, 10cm, 15cm

AEROCREPE™ Medium Cotton Crepe Bandage 2.5cm, 5cm, 7.5cm, 10cm, 15cm

AEROCREPE™ Elastic Cotton Crepe Bandage 2.5cm, 5cm, 7.5cm, 10cm, 15cm

It's recommended to have a selection of AEROCREPE<sup>™</sup> cotton bandage rolls to be prepared to address any relevant injuries.

#### Easy to open wrapper

All AEROCREPE $^{\text{TM}}$  cotton crepe bandages have an easy-open tab built into the wrapper of the bandage to allow for quick and easy use during an emergency, while keeping the contents safe and clean when unopened.

### **Contraindications and Precautions**

Any wound and gauze dressings underneath cotton compression bandages should be changed daily, or more often, according to the condition of the wound.

Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

### This Crepe Bandage product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Elastic Crepe Bandage 5cm