

# AEROCREPE LIGHT COTTON CREPE BANDAGE 7.5CM X 4M WRAP/12

ACL75





### **ADDITIONAL INFORMATION**

Weight 0.38 kg

**Dimensions**  $30 \times 5 \times 17 \text{ cm}$ 

Carton Qty 30

Pricing Unit Wrap/12

Size  $7.5 \text{cm} \times 4 \text{M}$ 

Min Order Qty 1

 ARTG ID
 236809

 GMDN Code
 10274

Medical Device Class 1

**UNSPSC** 42311500

**Tax Schedule** Taxable

Item Length 400

Item Width 7.5

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

Is Dangerous Good No

### **PRODUCT DESCRIPTION**

The AEROCREPE<sup>™</sup> Crepe Bandage is a high-quality cotton crepe-style bandage typically recommended for light to medium injury support requirements.

The crepe bandage comes in three different types (light, medium and elastic) to cater to the retention and protection of wound dressings, as well as provide assistance in standard immobilisation and compression bandage techniques.

## Protect gauze dressings with standard conforming bandages

The  $AEROCREPE^{TM}$  Light Cotton Crepe Bandage is of a standard conforming nature, with added durability that's suitable for wound dressing retention, as well as light injury support requirements.

These standard conforming bandages can be used to protect gauze bandages or gauze wraps placed around a wound and can easily



# AEROCREPE LIGHT COTTON CREPE BANDAGE 7.5CM X 4M WRAP/12

ACL75

adapt to support and strengthen the area of the injury.

Remember - protecting gauze dressings is paramount to avoiding infection in the wound.

### Comprehensive size ranges for all needs

Each AEROCREPE<sup>™</sup> cotton bandage roll comes in a length of 4 metres, with a variety of widths to suit all wound and injury requirements:

AEROCREPE<sup>™</sup> Light Cotton Crepe Bandage 5cm, 7.5cm, 10cm, 15cm

AEROCREPE™ Medium Cotton Crepe Bandage 2.5cm, 5cm, 7.5cm, 10cm, 15cm

AEROCREPE<sup>™</sup> Elastic Cotton Crepe Bandage 2.5cm, 5cm, 7.5cm,  $\frac{10cm}{10cm}$ 

It's recommended to have a selection of AEROCREPE<sup>™</sup> cotton bandage rolls to be prepared to address any relevant injuries.

# Easy to open wrapper

All AEROCREPE $^{\text{TM}}$  cotton crepe bandages have an easy-open tab built into the wrapper of the bandage to allow for quick and easy use during an emergency, while keeping the contents safe and clean when unopened.

### **Contraindications and Precautions**

Any wound and gauze dressings underneath cotton compression bandages should be changed daily, or more often, according to the condition of the wound.

Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

## This Crepe Bandage product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Light Crepe Bandage 7.5cm