

## **AEROCREPE MEDIUM COTTON CREPE BANDAGE** 10CM X 4M WRAP/12

**ACM100** 





### **ADDITIONAL INFORMATION**

Weight 0.5 kg

**Dimensions**  $27 \times 20 \times 4.5 \text{ cm}$ 

**Carton Qty** 20

**Pricing Unit** Wrap/12

Size 10cm x 4M

Min Order Qty

ARTG ID 236809 **GMDN Code** 10274

**Medical Device Class** Class 1

UNSPSC 42311500

**Tax Schedule** Taxable 400

**Item Width** 10

**Item Length** 

Item Measurement Unit cm

**Is Medical Device** Yes

Is Serialized Item No

**Is Dangerous Good** No

### **PRODUCT DESCRIPTION**

The AEROCREPE<sup>™</sup> Crepe Bandage is a high-quality cotton crepe-style bandage typically recommended for light to medium injury support requirements.

The crepe bandage comes in three different types (light, medium and elastic) to cater to the retention and protection of wound dressings, as well as provide assistance in standard immobilisation and compression bandage techniques.

## Protect gauze dressings with standard conforming bandages

The AEROCREPE™ Medium Cotton Crepe Bandage is also suitable for wound dressing retention, however, it's more applicable for medium injury support requirements.

These standard conforming bandages can be used to protect gauze bandages or gauze wraps placed around a wound and can easily



# AEROCREPE MEDIUM COTTON CREPE BANDAGE 10CM X 4M WRAP/12

**ACM100** 

adapt to support and strengthen the area of the injury.

Remember - protecting gauze dressings is paramount to avoiding infection in the wound.

### Comprehensive size ranges for all needs

Each AEROCREPE<sup>™</sup> cotton bandage roll comes in a length of 4 metres, with a variety of widths to suit all wound and injury requirements:

AEROCREPE<sup>™</sup> Light Cotton Crepe Bandage

5cm, 7.5cm, 10cm, 15cm

AEROCREPE<sup>™</sup> Medium Cotton Crepe Bandage

2.5cm, 5cm, 7.5cm, 10cm, 15cm

AEROCREPE<sup>™</sup> Elastic Cotton Crepe Bandage

2.5cm, 5cm, 7.5cm, 10cm, 15cm

It's recommended to have a selection of AEROCREPE™ cotton bandage rolls to be prepared to address any relevant injuries.

#### Easy to open wrapper

All AEROCREPE $^{\text{TM}}$  cotton crepe bandages have an easy-open tab built into the wrapper of the bandage to allow for quick and easy use during an emergency, while keeping the contents safe and clean when unopened.

#### **Contraindications and Precautions**

Any wound and gauze dressings underneath cotton compression bandages should be changed daily, or more often, according to the condition of the wound.

Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

### This Crepe Bandage product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Medium Crepe Bandage 10cm