



## ADDITIONAL INFORMATION

<b>Weight</b>	0.4 kg
<b>Dimensions</b>	27 × 15 × 4.5 cm
<b>Carton Qty</b>	30
<b>Pricing Unit</b>	Wrap/12
<b>Size</b>	7.5cm x 4M
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	236809
<b>GMDN Code</b>	10274
<b>Medical Device Class</b>	Class 1
<b>UNSPSC</b>	42311500
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	400
<b>Item Width</b>	7.5
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Barcode (GTIN/EAN/UPC)</b>	9341394001053

## PRODUCT DESCRIPTION

The AEROCREPE™ Crepe Bandage is a high-quality cotton crepe-style bandage typically recommended for light to medium injury support requirements.

The crepe bandage comes in three different types (light, medium and elastic) to cater to the retention and protection of wound dressings, as well as provide assistance in standard immobilisation and compression bandage techniques.

### Protect gauze dressings with standard conforming bandages

The AEROCREPE™ *Medium Cotton Crepe Bandage* is also suitable for wound dressing retention, however, it's more applicable for medium injury support requirements.

These standard conforming bandages can be used to protect gauze bandages or gauze wraps placed around a wound and can easily

adapt to support and strengthen the area of the injury.

Remember - protecting gauze dressings is paramount to avoiding infection in the wound.

### **Comprehensive size ranges for all needs**

Each AEROCREPE™ cotton bandage roll comes in a length of 4 metres, with a variety of widths to suit all wound and injury requirements:

*AEROCREPE™ Light Cotton Crepe Bandage*

[5cm](#), [7.5cm](#), [10cm](#), [15cm](#)

*AEROCREPE™ Medium Cotton Crepe Bandage*

[2.5cm](#), [5cm](#), [7.5cm](#), [10cm](#), [15cm](#)

*AEROCREPE™ Elastic Cotton Crepe Bandage*

[2.5cm](#), [5cm](#), [7.5cm](#), [10cm](#), [15cm](#)

It's recommended to have a selection of AEROCREPE™ cotton bandage rolls to be prepared to address any relevant injuries.

### **Easy to open wrapper**

All AEROCREPE™ cotton crepe bandages have an easy-open tab built into the wrapper of the bandage to allow for quick and easy use during an emergency, while keeping the contents safe and clean when unopened.

### **Contraindications and Precautions**

Any wound and gauze dressings underneath cotton compression bandages should be changed daily, or more often, according to the condition of the wound.

Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This Crepe Bandage product is not a substitute for direct medical attention from a healthcare practitioner.**

Medium Crepe Bandage 7.5cm