

## AEROFILM PLUS WATERPROOF ISLAND DRESSING 9 X 10CM BOX/20

AFPW009





## **ADDITIONAL INFORMATION**

Weight 0.12 kg

**Dimensions**  $17 \times 3.5 \times 13 \text{ cm}$ 

Carton Qty 8

Pricing Unit Box/20

Size 9 x 10cm

Min Order Qty 1

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311527

**Tax Schedule** Non-Taxable

Item Length10Item Width9

**Item Measurement Unit** cm

Is Medical Device Yes
Is Serialized Item No

**Is Dangerous Good** No

Shelf Life (Months) 36

## **PRODUCT DESCRIPTION**

AeroFilm Plus Waterproof Island Dressing provides transparent, waterproof protection with a soft absorbent island pad for comfortable post-operative and delicate wound coverage. The clear film acts as a bacterial barrier while allowing moisture vapour to escape, helping the site stay protected in daily life and letting users visually check the wound without removing the dressing.

The absorbent SupaSorb pad draws light exudate away from the skin and the low-adherent contact layer minimises trauma at change time. This combination supports a cleaner wound environment and gentler, low-pain removal—ideal for cuts, grazes, sutures and small procedural sites.

SupaBond hypoallergenic adhesive delivers reliable hold over body contours and joints, resisting edge lift from friction, sweat and movement. The film is ultra-thin for a discreet, conformable fit that moves with the body, maintaining a secure seal without bulk.



## AEROFILM PLUS WATERPROOF ISLAND DRESSING 9 X 10CM BOX/20

AFPW009

A unique frame delivery system simplifies accurate placement. The frame helps apply the dressing crease-free and ensures firm edge contact for a consistent seal in clinical settings and at home.

- Waterproof transparency: clear film shields the site while enabling visual checks.
- **Comfort absorbency:** SupaSorb pad manages exudate to help reduce maceration.
- **Hypoallergenic hold:** SupaBond adhesive provides strong, skin-friendly adhesion.
- Low-adherent contact: minimises sticking for gentler dressing changes.
- Easy, precise placement: frame delivery supports crease-free application.

**Size: 9 × 10 cm** offers broader coverage for small post-operative sites and sutures while remaining comfortable and conformable on contours, with transparent waterproof protection and a soft absorbent island pad.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Film Dressing Pad