



### ADDITIONAL INFORMATION

<b>Weight</b>	0.12 kg
<b>Dimensions</b>	17 × 3.5 × 13 cm
<b>Carton Qty</b>	8
<b>Pricing Unit</b>	Box/20
<b>Size</b>	9 x 10cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC</b>	42311527
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	10
<b>Item Width</b>	9
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36

### PRODUCT DESCRIPTION

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb™ absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.



# AEROFILM PLUS WATERPROOF ISLAND DRESSING

## 9 X 10CM BOX/20

AFPW009

Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™  
Hypo-allergenic adhesive  
exclusively engineered,  
for low irritation and high adhesion

SuperSorb™ Wound  
Exudate Pad  
Super Absorption,  
Reducing Skin Maceration

Non-stick Barrier  
Protects Wound with a Non-Adherent Barrier,  
Minimising Trauma on Removal

Dressing Frame

Ensures Easy, Effective Application

### **SupaBond™ Hypo-allergenic Adhesive**

- Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.

### **Transparent AeroFilm**

- The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.

### **Simple Application**

- AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over



# AEROFILM PLUS WATERPROOF ISLAND DRESSING

## 9 X 10CM BOX/20

AFPW009

the wound site.

### **Contraindications / Precautions**

- AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

**This product is not a substitute for direct medical attention from a healthcare practitioner.**

Film Dressing Pad