

# AEROFILM PLUS WATERPROOF ISLAND DRESSING 9 X 10CM BOX/3

**AFPW0093** 





#### **ADDITIONAL INFORMATION**

Weight 0.532 kg

**Dimensions**  $1.5 \times 13 \times 18.3 \text{ cm}$ 

Carton Qty 96

Pricing Unit Box/3

**Size** 9 x 10cm

Min Order Qty 1

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311527

**Tax Schedule** Non-Taxable

Item Length10Item Width9

Item Measurement Unit cm

**Is Medical Device** Yes

**Is Serialized Item** No

**Is Dangerous Good** No

36

Shelf Life (Months)

#### **PRODUCT DESCRIPTION**

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb<sup> $\mathsf{TM}$ </sup> absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.



### AEROFILM PLUS WATERPROOF ISLAND DRESSING 9 X 10CM BOX/3

**AFPW0093** 

AeroFilm Ultra-thin, Waterproof & Transparent Bacterial Barrier
Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion
SuperSorb™ Wound Exudate Pad Super Absorption, Reducing Skin Maceration
Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal
Dressing Frame
Ensures Easy, Effective Application
SupaBond™ Hypo-allergenic Adhesive
<ul> <li>Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring th dressing stays in place over body contours.</li> </ul>
Transparent AeroFilm

• The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.



# AEROFILM PLUS WATERPROOF ISLAND DRESSING Healthcare

**AFPW0093** 

Simple Applicati	ion
the	roFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables e dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over e wound site.
Contraindication	ns / Precautions
∘ Ae	roFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
∘ As	with all adhesive products, apply and remove carefully from sensitive or fragile skin.
∘ See coi	ek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and nsult a physician.



## AEROFILM PLUS WATERPROOF ISLAND DRESSING 9 X 10CM BOX/3

**AFPW0093** 

This product is not a substitute for direct medical attention from a healthcare practit
---

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Film Dressing Pad