



### ADDITIONAL INFORMATION

<b>Weight</b>	0.15 kg
<b>Dimensions</b>	19.7 × 3.5 × 14 cm
<b>Carton Qty</b>	8
<b>Pricing Unit</b>	Box/20
<b>Size</b>	10cm x 12cm
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	159198
<b>GMDN Code</b>	44990
<b>Medical Device Class</b>	Class 1s
<b>UNSPSC Code</b>	42311527
<b>Tax Schedule</b>	Non-Taxable
<b>Item Length</b>	12
<b>Item Width</b>	10
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No
<b>Shelf Life (Months)</b>	36
<b>Warranty Period (Months)</b>	12



### PRODUCT DESCRIPTION

AeroFilm Plus Waterproof Island Dressing pairs a transparent, waterproof film with a soft, absorbent island pad to protect delicate and

post-operative sites through everyday activity. The clear polyurethane-style film forms a barrier to water and external contaminants while allowing moisture vapour to pass, helping the skin feel fresher under the dressing. Because the site remains visible through the film, quick visual checks can be made without disturbing the seal. Each dressing is sterile and single-use, making it suitable for clinical care, workplace first aid stations, and active households.

At the centre, the SupaSorb absorbent pad manages light exudate to help reduce maceration. A low-adherent contact layer minimises sticking to fragile tissue so removal is gentler and less disruptive. This balanced design supports a cleaner environment around cuts, grazes, sutures and small procedural sites, while the low profile sits comfortably under clothing and gloves.

SupaBond hypoallergenic adhesive provides dependable hold over contours and joints. It resists edge lift from friction, perspiration and motion to help the dressing stay in place during showering, work shifts and daily routines. The ultra-thin, conformable film moves smoothly with the body for a discreet finish that supports confidence during wear. The dressing is latex free and suitable for sensitive skin.

An integrated frame delivery system guides accurate placement. The frame helps apply the dressing crease-free and ensures firm edge contact for a consistent, reliable seal. This simple, intuitive application is valued in clinics for speed and repeatability and at home for ease of use and consistent results. From minor procedures to routine first aid, AeroFilm Plus delivers transparent, waterproof protection with comfortable absorbency.

- **Waterproof transparency:** clear film shields the site and enables visual checks.
- **Comfort absorbency:** SupaSorb pad manages exudate to help reduce maceration.
- **Hypoallergenic strong hold:** SupaBond adhesive provides reliable adhesion on contours.
- **Low-adherent contact:** minimises sticking for gentler, low-pain removal.
- **Easy, precise placement:** frame delivery supports smooth, crease-free application.

The **Size: 10 x 12 cm** option offers broader coverage for larger post-operative sites or dressings placed over high-movement areas such as the lower back or hip, while remaining thin and conformable for comfortable wear under clothing.

Supplied in a **box of 20 dressings**, ideal for clinics, workplace cabinets and first aid kits where regular use requires convenient on hand stock. This pack size supports efficient replenishment and helps keep dependable waterproof film dressings readily available.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.