

AEROFILM PLUS WATERPROOF ISLAND DRESSING 15 X 20CM BOX/20

AFPW015





ADDITIONAL INFORMATION

Weight 0.33 kg

Dimensions $30 \times 3.5 \times 18.5 \text{ cm}$

Carton Qty 8

Pricing Unit Box/20

Size 15cm x 20cm

Min Order Qty 1

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311527

Tax Schedule Non-Taxable

Item Length20Item Width15Item Measurement Unitcm

Is Medical Device Yes

Is Serialized Item No

Is Dangerous Good No

Shelf Life (Months) 36

PRODUCT DESCRIPTION

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb™ absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.



AEROFILM PLUS WATERPROOF ISLAND DRESSING 15 X 20CM BOX/20

AFPW015

AeroFilm Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad Super Absorption, Reducing Skin Maceration

Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal

Dressing Frame

Ensures Easy, Effective Application

SupaBond™ Hypo-allergenic Adhesive

• Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.

Transparent AeroFilm

• The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.

Simple Application

• AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site.

Contraindications / Precautions

- AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Film Dressing Pad



AEROFILM PLUS WATERPROOF ISLAND DRESSING

AFPW015