

# AEROFILM PLUS WATERPROOF ISLAND DRESSING 4 X 5CM BOX/3

**AFPW0043** 





#### **ADDITIONAL INFORMATION**

Weight 0.02 kg

**Dimensions**  $1.5 \times 9.5 \times 15 \text{ cm}$ 

Carton Qty 96

Pricing Unit Box/3

Size 4cm x 5cm

Min Order Qty 1

**ARTG ID** 159198 **GMDN Code** 44990

Medical Device Class 1s

**UNSPSC** 42311527

**Tax Schedule** Non-Taxable

No

Item Length5Item Width4

Item Measurement Unit cm

**Is Medical Device** Yes

Is Serialized Item

**Is Dangerous Good** No

Shelf Life (Months) 36

#### **PRODUCT DESCRIPTION**

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb<sup> $\mathsf{TM}$ </sup> absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.



## AEROFILM PLUS WATERPROOF ISLAND DRESSING 4 X 5CM BOX/3

**AFPW0043** 

AeroFilm Ultra-thin, Waterproof & Transparent Bacterial Barrier

Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion

SuperSorb™ Wound Exudate Pad Super Absorption, Reducing Skin Maceration

Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal

**Dressing Frame** 

Ensures Easy, Effective Application

#### SupaBond™ Hypo-allergenic Adhesive

• Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.

#### **Transparent AeroFilm**

• The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without removing the dressing. This allows for maximum usage of the dressing before removal.

#### **Simple Application**

• AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site.

### **Contraindications / Precautions**

- AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

## This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Film Dressing Pad



# AEROFILM PLUS WATERPROOF ISLAND DRESSING AVECUA DOM:

**AFPW0043**