

AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/20

AFPW00620





ADDITIONAL INFORMATION

Weight 0.77 kg

Dimensions $15 \times 3.5 \times 10 \text{ cm}$

Carton Qty 36

Pricing Unit Box/20

Size 6cm x 7cm

Min Order Qty 1

ARTG ID 159198 **GMDN Code** 44990

Medical Device Class 1s

UNSPSC 42311527

Tax Schedule Non-Taxable

Item Length7Item Width6

Item Measurement Unit cm

Is Medical Device Yes
Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

This post-operative dressing combines the advantages of AEROFILM with a SupaSorb^{\mathbf{m}} absorbent and low adherent dressing pad to provide a superior dressing for post operative and delicate wound sites.

AeroFilm



AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/20

AFPW00620

Ultra-thin, Waterproof & Transparent Bacterial Barrier
Supabond™ Hypo-allergenic adhesive exclusively engineered, for low irritation and high adhesion
SuperSorb™ Wound Exudate Pad Super Absorption, Reducing Skin Maceration
Non-stick Barrier Protects Wound with a Non-Adherent Barrier, Minimising Trauma on Removal
Dressing Frame
Ensures Easy, Effective Application
SupaBond™ Hypo-allergenic Adhesive
 Exclusively engineered for unrivalled adhesion, this hypoallergenic adhesive has a low irritation risk while ensuring the dressing stays in place over body contours.
Transparent AeroFilm

o The transparent AeroFilm allows the user to monitor the level of fluid and check for clinical signs of infection without

removing the dressing. This allows for maximum usage of the dressing before removal.



AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/20

AFPW00620

Simple Application
 AeroFilm has an integrated frame applicator that provides a simple and effective application system. The frame enables the dressing to be applied crease free and ensures adhesion. AeroFilm transparency enables correct positioning over the wound site.
Contraindications / Precautions
o AeroFilm should not be used on infected wounds unless under the supervision of a healthcare practitioner.
○ As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
 Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.



AEROFILM PLUS WATERPROOF ISLAND DRESSING 6 X 7CM BOX/20

AFPW00620

This product is not a substitute for direct medical attention from a healthcare practitioner.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

CONTENTS LISTS

Code	Description	Qty
AFPW006	AEROFILM PLUS Waterproof Island Dressing 6 x 7cm	1

FILM DRESSING PAD